


MoliNea® & Strampelpeter®

HARTMANN




Fitting guide: Rectangular Pads and Underpads


MoliNea®, MoliNea® Couche & Strampelpeter®



Use **MoliNea®** rectangular pads in conjunction with **MoliCare® Premium Fixpants®** or close fitting underwear.




Place inside **MoliCare® Premium Fixpants®** or underwear coloured side down.



Stick the adhesive strip into the underwear to secure in place. (applicable only for rectangular pads with adhesive strip)

MoliNea® Underpads



Place on bed or other surface coloured side down.

MoliNea®

- Waterproof backing • Ideal for faecal incontinence



For very light to moderate bladder weakness.



Code	Absorbency Level	Working Capacity	Pack Quantity	Case Quantity
160 890	● 1 drop	140ml	50	6
160 892	● 1 drop	207ml	50	4
160 850	● 1 drop	387ml	28	4
160 885	● 1 drop	347ml	30	4

MoliNea® Couche

- Waterproof backing • Designed to be worn within underwear



For very light to moderate bladder weakness.



Code	Absorbency Level	Working Capacity	Pack Quantity	Case Quantity
148 510	● 1 drop	117ml	56	4
148 520	● 1 drop	161ml	56	4
148 530	● 1 drop	241ml	56	4

Strampelpeter®

- Insert pad ideally designed to be used inside pouch pants



For very light bladder weakness.




Code	Absorbency Level	Working Capacity	Pack Quantity	Case Quantity
148 610	● 1 drop	147ml	56	4
148 620	● 1 drop	230ml	56	4

MoliNea® Underpad

- Stay-dry, non woven top layer for patient comfort • Waterproof backing



For light bladder weakness.



Code	Absorbency Level	Working Capacity	Dimensions (cm)	Pack Quantity	Case Quantity
161 221	● 2 drops	101ml	40 x 60	30	6
161 320	● 2 drops	140ml	60 x 60	30	4
161 520	● 2 drops	221ml	60 x 90	30	4
161 100	● 2 drops	169ml	40 x 60	150	1
161 182	● 2 drops	288ml	60 x 60	25	4
161 600	● 2 drops	400ml	60 x 90	100	1
161 804	● 2 drops	316ml	90 x 60 (with wings)	20	4
163 300	● 3 drops	400ml	60 x 60	100	1
161 194	● 3 drops	625ml	60 x 90	25	4

To order additional products, go to www.hartmann-direct.co.uk or call 0800 028 9490